R Dose rates in µSv/h from the pixel detector after 10 years of running and 15 days of cooling. →									
1 m - •5	•5	•4	•3	•2	•2	•1	•1	•1	•1
•8	•8	•6	•4	•3	•2	•1	•1	•1	•1
0.5 m +15	•15	•11	•6	•3	•2	•2	•2	•2	•2
•20	•21	•15	•7	•4	•2	•2	•2	•3	•2
•29	•33	•22	•9	•4	•2	•2	•2	•5	•2
•40	•51	•32 PP	•12 201	●4	•2	•2	•2	•9	•3
70		•37	•13	•5	•3	•2	•3 120	•19	•3
		•35	•14	•5	•3	•3	•3 19		•3
		•37	•14	•7	•5	•4	•5 460		-
					Beamp				7
		1			1	*		<u>_</u>]-]	4 m